


FULL DENTURES

Dentures Can Look Artificial or They Can Look Simply Natural!



Old Denture



New Denture



Old Denture



New Denture




Old Denture



New Denture

Ask your dentist how Simply Natural Dentures can improve your smile.



**Simply
Natural
Dentures®**

Glidewell Laboratories

©2003 Glidewell Laboratories

Photos courtesy of Geneva Dental Inc., Beverly Hills, CA

01-1167-0603




Simply Natural Dentures®

Your dentist is among a select few who use the Simply Natural Dentures system for more natural looking full dentures. They resist staining, are stronger, last longer, and fit more comfortably than typical dentures. Here are the easy steps required to make your new Simply Natural Dentures.

At your **first appointment** your dentist will ask you some specific questions about your present dentures, and how you want your new dentures to look. Impressions are taken and the models sent to the lab for fabrication of your esthetic guide.



Esthetic Guide



Old



New

Lip Support

During the **second appointment** your doctor uses a unique esthetic guide to gather your personal smile characteristics. Your dentist positions it in your mouth, sets your bite, and readjusts the two front teeth if necessary. This helps guide the dental technician in creating your new smile. Special consideration is taken to ensure proper lip support and a pleasing smile line for your facial contour.



Kenson Teeth



Premium Teeth

At this same appointment the specific shade and shape of the denture teeth are selected. For a more natural looking denture, hand sculpted European denture teeth can be selected at an additional cost. These teeth give the most natural appearance and are virtually indistinguishable from real teeth.



At the laboratory, the technician uses your esthetic guide as a foundation to carefully position your new teeth. Your new dentures are contoured in a manner which regains the appearance of natural healthy tissue.

At the **third appointment** a wax copy of your new denture is tried in to confirm bite and tooth position. Any changes to your smile or bite can be easily made at this point. When you're fully satisfied your denture is sent to the laboratory for final processing.

In the laboratory, we process your new dentures using a system that ensures a comfortable fit. A special bond enhancer is placed on the teeth to secure them permanently into the denture. After processing, your dentures are shaped to look like natural healthy tissue making them comfortable and easy to clean.



Before



After

At your **delivery appointment** you will receive your new dentures. You'll find your new Simply Natural Dentures to be the most lifelike and natural dentures you've ever had.

"My new dentures give me more confidence in everything I do."

